

Youcare Ostomy System




Feature 1

Assurance for Skin Friendly

Skin-Friendly

Skin Barrier of Natural Materials




The main ingredients of the skin barrier are the natural materials that rarely irritate the skin such as “karaya gum” and “citrus pectin”. Karaya gum keeps the surface of the faceplate weak acid.

Feature 2


Easy and Simple Handling

User-Friendly

Clip-less Ostomy Appliances



Just turn up the plate along with the blue fold up line and secure the end with Velcro closures attached. Opening and closing are easy. The keep-shaping plate makes discharge and wiping easy.




By folding the keep-shaping plate in V shape, open-end can be kept open.

Feature 3

Comfortable Application

Scientific Analysis of Body Movement

Diamond Shape



The oval shape flattens the fine wrinkles over the abdomen.

The triangle shape achieves excellent fitting to the inguinal region.

Anatomically designed Diamond Shape effectively prevents detachment and leakage by increased flexibility following the curved groin and the abdominal movement. Because of the compact design, he/she hardly notices uncomfortable feeling. It offers higher levels of durability and achieves excellent contact.

Youcare Ostomy System

Digestive System Ostomy Appliances

1-Piece Ostomy System



Open-ended Clip
Flat

Youcare-D

2 Days Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
Flesh	16761	20mm ϕ	14~59mm (free-cut)	275mm X 150mm	30pcs
Without Filter	16771	20mm ϕ	14~59mm (free-cut)		
20	16781	20mm ϕ	14~59mm (free-cut)		
25	16782	25mm ϕ	19~24mm		
30	16783	30mm ϕ	24~29mm		
35	16784	35mm ϕ	29~34mm		
40	16785	40mm ϕ	34~39mm		
45	16786	45mm ϕ	39~44mm		
50	16787	50mm ϕ	44~49mm		




Open-ended Clip
Convex

Youcare-Dc

2 Days Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
22	16811	22mm ϕ	18~21mm	275mm X 150mm	10pcs
25	16812	25mm ϕ	21~24mm		
28	16813	28mm ϕ	24~27mm		
32	16814	32mm ϕ	27~31mm		
36	16815	36mm ϕ	31~35mm		
40	16816	40mm ϕ	35~39mm		
45	16817	45mm ϕ	39~44mm		
50	16818	50mm ϕ	44~49mm		



Open-ended Clipless
Flat

Youcare-TD

2 Days Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
Flesh	17441	20mm ϕ	14~59mm (free-cut)	275mm X 150mm	10pcs
Without Filter	17451	20mm ϕ	14~59mm (free-cut)		
20	17461	20mm ϕ	14~59mm (free-cut)		
25	17462	25mm ϕ	19~24mm		
30	17463	30mm ϕ	24~29mm		
35	17464	35mm ϕ	29~34mm		
40	17465	40mm ϕ	34~39mm		
45	17466	45mm ϕ	39~44mm		
50	17467	50mm ϕ	44~49mm		



Open-ended Clipless
Convex

Youcare-TDc

2 Days Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
22	17471	22mm ϕ	18~21mm	275mm X 150mm	10pcs
25	17472	25mm ϕ	21~24mm		
28	17473	28mm ϕ	24~27mm		
32	17474	32mm ϕ	27~31mm		
36	17475	36mm ϕ	31~35mm		
40	17476	40mm ϕ	35~39mm		
45	17477	45mm ϕ	39~44mm		
50	17478	50mm ϕ	44~49mm		



Closed
Flat

Youcare-C

1 Day Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
Flesh	16791	20mm ϕ	14~59mm (free-cut)	200mm X 150mm	30pcs
20	16801	20mm ϕ	14~59mm (free-cut)		
25	16802	25mm ϕ	19~24mm		
30	16803	30mm ϕ	24~29mm		
35	16804	35mm ϕ	29~34mm		
40	16805	40mm ϕ	34~39mm		
45	16806	45mm ϕ	39~44mm		
50	16807	50mm ϕ	44~49mm		




Closed
Convex

Youcare-Cc

1 Day Exchange

Type	Product Code No.	Specification		Bag Size (Length X Width)	Contents Per Box
		Starter Hole Size	Max Available Size		
22	16821	22mm ϕ	18~21mm	200mm X 150mm	10pcs
25	16822	25mm ϕ	21~24mm		
28	16823	28mm ϕ	24~27mm		
32	16824	32mm ϕ	27~31mm		
36	16825	36mm ϕ	31~35mm		
40	16826	40mm ϕ	35~39mm		
45	16827	45mm ϕ	39~44mm		
50	16828	50mm ϕ	44~49mm		

2-Piece Ostomy System



Flat

Youcare 2-F

3 to 5 Days Exchange

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M22	16591	22mm ϕ	16~34mm (free-cut)	97mmX105mm	M	5pcs
M25	16592	25mm ϕ	21~24mm	97mmX105mm	M	
M28	16593	28mm ϕ	24~27mm	97mmX105mm	M	
L32	16594	32mm ϕ	27~31mm	108mmX117mm	L	
L36	16595	36mm ϕ	31~35mm	108mmX117mm	L	
L40	16596	40mm ϕ	35~39mm	108mmX117mm	L	
L45	16597	10mm ϕ	4~44mm (free-cut)	108mmX117mm	L	
65	16598	10mm ϕ	4~64mm (free-cut)	129mmX140mm	65	



Convex

Youcare 2-Fc

3 to 5 Days Exchange

Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size	Applicable Bag	
M22	16601	22mm ϕ	18~21mm	110mm ϕ	M	5pcs
M25	16602	25mm ϕ	21~24mm	110mm ϕ	M	
M28	16603	28mm ϕ	24~27mm	110mm ϕ	M	
L32	16604	32mm ϕ	27~31mm	120mm ϕ	L	
L36	16605	36mm ϕ	31~35mm	120mm ϕ	L	
L40	16606	40mm ϕ	35~39mm	120mm ϕ	L	

* C : closed c : convex D : drainable F : flange f : filter TD : tape drainable U : urostomy

2-Piece Ostomy System



Open-ended Clip

Youcare 2-D

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
Clear	M 16621	300mmX150mm		10pcs
	L 16622	300mmX150mm		
	65 16623	300mmX150mm		
Flesh	M 17751	300mmX150mm		
	L 17752	300mmX150mm		



Open-ended Clip

Youcare 2-Df

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
Clear	M 16611	300mmX150mm		10pcs
	L 16612	300mmX150mm		
	65 16613	300mmX150mm		
Flesh	M 17741	300mmX150mm		
	L 17742	300mmX150mm		



Open-ended Clipless

Youcare 2-TD

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
Clear	M 17781	300mmX150mm		10pcs
	L 17782	300mmX150mm		
	65 17783	300mmX150mm		
Flesh	M 17791	300mmX150mm		
	L 17792	300mmX150mm		



Open-ended Clipless

Youcare 2-TDf

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
Clear	M 17761	300mmX150mm		10pcs
	L 17762	300mmX150mm		
	65 17763	300mmX150mm		
Flesh	M 17771	300mmX150mm		
	L 17772	300mmX150mm		



Closed

Youcare 2-Cf

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
Clear	M 16631	210mmX150mm		10pcs
	L 16632	210mmX150mm		
Flesh	M 17801	210mmX150mm		
	L 17802	210mmX150mm		



Closed

Youcare 2-M

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
M	16671	120mmX130mm		10pcs
L	16672	120mmX130mm		10pcs

Urology System Ostomy Appliances

1-Piece Ostomy System




Outlet for Urine
Flat

Youcare-U

2 to 3 Days Exchange

Type	Product Code No.	Starter Hole Size	Specification		Bag Size (Length X Width)	Contents Per Box
			Max Available Size			
10	15361	10mm ϕ	4~59mm (free-cut)		260mm X 150mm	10pcs
20	15362	20mm ϕ	14~19mm			
25	15363	25mm ϕ	19~24mm			
30	15364	30mm ϕ	24~29mm			
35	15365	35mm ϕ	29~34mm			
40	15366	40mm ϕ	34~39mm			



Outlet for Urine
Convex


Youcare-Uc

2 to 3 Days Exchange

Type	Product Code No.	Starter Hole Size	Specification		Bag Size (Length X Width)	Contents Per Box
			Max Available Size			
13	16401	13mm ϕ	9~12mm		260mm X 150mm	10pcs
16	16402	16mm ϕ	12~15mm			
19	16403	19mm ϕ	15~18mm			
22	16404	22mm ϕ	18~21mm			
25	16405	25mm ϕ	21~24mm			
28	16406	28mm ϕ	24~27mm			
32	16407	32mm ϕ	27~31mm			

Urology System Ostomy Appliances

2-Piece Ostomy System




Flat

Youcare 2-Fu

4 to 7 Days Exchange

Type	Product Code No.	Specification				Applicable Bag	Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size			
M13	16641	13mm ϕ	7~34mm (free-cut)	97mmX105mm	M	5pcs	
M16	16642	16mm ϕ	12~15mm	97mmX105mm	M		
M19	16643	19mm ϕ	15~18mm	97mmX105mm	M		
M22	16644	22mm ϕ	18~21mm	97mmX105mm	M		
M25	16645	25mm ϕ	21~24mm	97mmX105mm	M		
M28	16646	28mm ϕ	24~27mm	97mmX105mm	M		
L32	16647	32mm ϕ	27~44mm (free-cut)	108mmX117mm	L		



Convex

Youcare 2-Fuc

4 to 7 Days Exchange

Type	Product Code No.	Specification				Applicable Bag	Contents Per Box
		Starter Hole Size	Max Available Size	Face Plate Size			
M13	16651	13mm ϕ	9~12mm	110mm ϕ	M	5pcs	
M16	16652	16mm ϕ	12~15mm	110mm ϕ	M		
M19	16653	19mm ϕ	15~18mm	110mm ϕ	M		
M22	16654	22mm ϕ	18~21mm	110mm ϕ	M		
M25	16655	25mm ϕ	21~24mm	110mm ϕ	M		
M28	16656	28mm ϕ	24~27mm	110mm ϕ	M		
L32	16657	32mm ϕ	27~31mm	120mm ϕ	L		



Outlet for Urine

Youcare 2-U

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
M	16661	260mmX150mm		10pcs
L	16662	260mmX150mm		10pcs

Accessories



Youcare Trunk Belt

Type	Product Code No.	Specification		Contents Per Box
		Waist Size	Applicable Bag	
—	15091	60~110mm	Youcare 2-TDf Youcare 2-TD Youcare 2-Df Youcare 2-D Youcare 2-Cf Youcare 2-U Youcare 2-IS Youcare 2-M Youcare 2-BC Youcare-TDc Youcare-Dc Youcare-Cc Youcare-Uc Ilefine D-Cap Ilefine TD Ilefine D	1Pce




Youcare Trunk Belt With Face Plate

Type	Product Code No.	Specification		Contents Per Box
		Waist Size	Applicable Bag	
—	15081	60~110mm (diameter of hole 73mm ϕ)	Youcare-TD Youcare-D Youcare-C Youcare-U Pospack-B Pospack-Light	1set



Youcare 2-BC

Type	Product Code No.	Specification (Height of Dome)		Contents Per Box
M	16681	20mm		2pcs
L	16682	20mm		2pcs
65	16683	20mm		2pcs



Youcare 2-IS

Type	Product Code No.	Specification		Contents Per Box
		Bag Size (Length X Width)		
M	16691	900mmX100mm		5pcs
L	16692	900mmX150mm		5pcs



ALCARE Co., Ltd.

International Sales Div.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan

TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp

www.alcare.co.jp/e