

# Ostomy Accessories

ALCARE provides the skin barrier containing the natural raw materials with excellent skin protecting and a large variety of pouch accessories so that the ostomates can enjoy more comfortable daily life.

## Skin Barrier




### Procare Wafer

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- It can be freely cut to the shape of the hollow or the wrinkled area, or it can be used in layers if necessary.

Type	Product Code No.	Specification (Length X Width)	Contents Per Box
10	12581	100mmX100mm	10pcs
20	12582	200mmX195mm	5pcs



### Procare Soft Wafer-Ring

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- A moderately soft wafer that can be molded with the hands (no scissors needed).
- A ring-shaped wafer that is easy to apply to the peristomal skin area.
- A double-sided adhesive wafer that adheres well to the face plate of various appliances.

Type	Product Code No.	Specification (Diameter of Hole)	Contents Per Box
Ring	20	20mmφ	8pcs
	30	30mmφ	8pcs



### Procare Soft Wafer-Stick

Wafer Skin Barrier (Karaya Gum + Citrus Pectin)

- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- An extremely soft wafer that can be freely molded.
- A double-sided adhesive wafer that adheres well to the face plate of various appliances.
- Ideal for application to uneven hard-to-apply surfaces such as the area around the fistula.

Type	Product Code No.	Specification (Length X Width)	Contents Per Box
Stick	15601	15mmX100mm	10pcs



### Procure MF Pate

Paste Skin Barrier (Karaya Gum + Citrus Pectin)

- It fills the gap between the stoma and face plate, protecting the skin from discharge.
- An alcohol-free paste which will not make the skin sting but will remain soft.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12601	100g	1tube

# Ostomy Accessories



## Procure Paste

Paste Skin Barrier (Karaya Gum + Citrus Pectin)

- By applying this paste on the inside of the face plate, the adhesion of appliances can be improved.
- A laminate tube keeps the paste moist, making it easy to squeeze the tube down to the last.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	14431	50g	1tube



## Procure Powder

Powder Skin Barrier (Karaya Gum + Citrus Pectin)

- This powder fills the gap between the stoma and the face plate, protecting the skin from discharge.
- It is easy to apply, even to a small gap. The powder offers excellent water absorption.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12611	50g	1bottle

## Stoma Care Accessory



## Procure Remover

Ostomy Appliance Remover

- It reduces the peeling pain when removing the face plate of ostomy appliance from the skin.
- It can be effectively used when an appliance needs immediate removal.
- It is individually packed in a foil pack.

Type	Product Code No.	Specification (Fluid Measure)	Contents Per Box
—	13871	2mL	50pcs



## Deofine Powder

Powdery Deodorant for Ostomates

- Plant ingredient inhibits the generation and diffusion of odors.
- Refreshing mint ingredient makes stoma care more comfortable.
- Stick-type individual packaging can provide portability.

Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	17891	1.5g / pack	30pks



**ALCARE Co., Ltd.**

International Sales Div.

19F, Arcacentral Bldg., 1-2-1 Kinshi, Sumida-ku, Tokyo 130-0013, Japan

TEL.+81-3-5611-7819 FAX.+81-3-5611-7827 MAIL. info@alcare.co.jp

www.alcare.co.jp/e