


Youcare-TD

1-Piece Ostomy System

Features

- Elastic and durable thin skin barrier is used for faceplate.
- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool (Without-filter type is also available).
- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L × W)	
Flesh	17441	20 mmΦ	14-59 mm (Free-cut)	275 mm × 150 mm	10 pieces
W/O Filter	17451	20 mmΦ	14-59 mm (Free-cut)	275 mm × 150 mm	10 pieces
20	17461	20 mmΦ	14-59 mm (Free-cut)	275 mm × 150 mm	10 pieces
25	17462	25 mmΦ	19-24 mm	275 mm × 150 mm	10 pieces
30	17463	30 mmΦ	24-29 mm	275 mm × 150 mm	10 pieces
35	17464	35 mmΦ	29-34 mm	275 mm × 150 mm	10 pieces
40	17465	40 mmΦ	34-39 mm	275 mm × 150 mm	10 pieces
45	17466	45 mmΦ	39-44 mm	275 mm × 150 mm	10 pieces
50	17467	50 mmΦ	44-49 mm	275 mm × 150 mm	10 pieces