


Youcare-TDc

1-Piece Ostomy System

Features

- Elastic and durable thin skin barrier is used for faceplate.
- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.
- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L × W)	
22	17471	22 mmΦ	18-21 mm	275 mm × 150 mm	10 pieces
25	17472	25 mmΦ	21-24 mm	275 mm × 150 mm	10 pieces
28	17473	28 mmΦ	24-27 mm	275 mm × 150 mm	10 pieces
32	17474	32 mmΦ	27-31 mm	275 mm × 150 mm	10 pieces
36	17475	36 mmΦ	31-35 mm	275 mm × 150 mm	10 pieces
40	17476	40 mmΦ	35-39 mm	275 mm × 150 mm	10 pieces
45	17477	45 mmΦ	39-44 mm	275 mm × 150 mm	10 pieces
50	17478	50 mmΦ	44-49 mm	275 mm × 150 mm	10 pieces