

Youcare-D

1-Piece Ostomy System

Features

- Elastic and durable thin skin barrier is used for faceplate.
- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool (Without-filter type is also available).



| Type | Product Code No. | Specification | | | Contents Per Box |
|------------|------------------|-------------------|------------------------|------------------|------------------|
| | | Starter Hole Size | Max Available Size | Bag Size (L × W) | |
| Flesh | 16761 | 20 mmΦ | 14-59 mm (Free-cut) | 275 mm × 150 mm | 30 pieces |
| W/O Filter | 16771 | 20 mmΦ | 14-59 mm (Free-cut) | 275 mm × 150 mm | 30 pieces |
| 20 | 16781 | 20 mmΦ | 14-59 mm (Free-cut) | 275 mm × 150 mm | 30 pieces |
| 25 | 16782 | 25 mmΦ | 19-24 mm | 275 mm × 150 mm | 30 pieces |
| 30 | 16783 | 30 mmΦ | 24-29 mm | 275 mm × 150 mm | 30 pieces |
| 35 | 16784 | 35 mmΦ | 29-34 mm | 275 mm × 150 mm | 30 pieces |
| 40 | 16785 | 40 mmΦ | 34-39 mm | 275 mm × 150 mm | 30 pieces |
| 45 | 16786 | 45 mmΦ | 39-44 mm | 275 mm × 150 mm | 30 pieces |
| 50 | 16787 | 50 mmΦ | 44-49 mm | 275 mm × 150 mm | 30 pieces |