

Youcare-C

1-Piece Ostomy System

Features

- Elastic ultra-thin skin barrier is used for faceplate.
- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L × W)	
Flesh	16791	20 mmΦ	14-59 mm (Free-cut)	200 mm × 150 mm	30 pieces
20	16801	20 mmΦ	14-59 mm (Free-cut)	200 mm × 150 mm	30 pieces
25	16802	25 mmΦ	19-24 mm	200 mm × 150 mm	30 pieces
30	16803	30 mmΦ	24-29 mm	200 mm × 150 mm	30 pieces
35	16804	35 mmΦ	29-34 mm	200 mm × 150 mm	30 pieces
40	16805	40 mmΦ	34-39 mm	200 mm × 150 mm	30 pieces
45	16806	45 mmΦ	39-44 mm	200 mm × 150 mm	30 pieces
50	16807	50 mmΦ	44-49 mm	200 mm × 150 mm	30 pieces