


Youcare 2-F

2-Piece Ostomy System

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Thinner outer edge of faceplate assures close contact with body and reduces uncomfortableness while attached.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M22	16591	22 mmΦ	16-34 mm (Free-cut)	97 mm × 105 mm	M	5 pieces
M25	16592	25 mmΦ	21-24 mm	97 mm × 105 mm	M	5 pieces
M28	16593	28 mmΦ	24-27 mm	97 mm × 105 mm	M	5 pieces
L32	16594	32 mmΦ	27-31 mm	108 mm × 117 mm	L	5 pieces
L36	16595	36 mmΦ	31-35 mm	108 mm × 117 mm	L	5 pieces
L40	16596	40 mmΦ	35-39 mm	108 mm × 117 mm	L	5 pieces
L45	16597	10 mmΦ	4-44 mm (Free-cut)	108 mm × 117 mm	L	5 pieces
65	16598	10 mmΦ	4-64 mm (Free-cut)	129 mm × 140 mm	65	5 pieces