

# Youcare 2-TDf

## 2-Piece Ostomy System

### Features

- Highly functional deodorizing filter releases gas smoothly while preventing leakage of odor or stool.
- Clip-less model provides you with simple operation and reduction of uncomfortableness while attached.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type		Product Code No.	Specification	Contents Per Box
			Bag Size (L x W)	
Transparent	M	17761	300 mm × 150 mm	10 pieces
	L	17762	300 mm × 150 mm	10 pieces
	65	17763	300 mm × 150 mm	10 pieces
Flesh	M	17771	300 mm × 150 mm	10 pieces
	L	17772	300 mm × 150 mm	10 pieces