

Youcare 2-D

2-Piece Ostomy System

Features

- Stoma bag without a deodorizing filter.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type		Product Code No.	Specification	Contents Per Box
			Bag Size (L × W)	
Transparent	M	16621	300 mm × 150 mm	10 pieces
	L	16622	300 mm × 150 mm	10 pieces
	65	16623	300 mm × 150 mm	10 pieces
Flesh	M	17751	300 mm × 150 mm	10 pieces
	L	17752	300 mm × 150 mm	10 pieces