

Youcare-Uc

1-Piece Ostomy System

Features

- Elastic and durable thin skin barrier, which is specially developed for the urostomy, is used for faceplate.
- Youcare Trunk Belt can improve close contact between faceplate and skin.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (L × W)	
13	16401	13 mmΦ	9-12 mm	260 mm × 150 mm	10 pieces
16	16402	16 mmΦ	12-15 mm	260 mm × 150 mm	10 pieces
19	16403	19 mmΦ	15-18 mm	260 mm × 150 mm	10 pieces
22	16404	22 mmΦ	18-21 mm	260 mm × 150 mm	10 pieces
25	16405	25 mmΦ	21-24 mm	260 mm × 150 mm	10 pieces
28	16406	28 mmΦ	24-27 mm	260 mm × 150 mm	10 pieces
32	16407	32 mmΦ	27-31 mm	260 mm × 150 mm	10 pieces