


Youcare 2-Fu

2-Piece Ostomy System

Features

- Recommended for people who favor without tapes as skin barrier is used on entire surface.
- Thinner outer edge of faceplate assures close contact with body and reduces uncomfortableness while attached.



Type	Product Code No.	Specification				Contents Per Box
		Starter Hole Size	Max Available Size	Faceplate Size	Applicable Bag	
M13	16641	13 mmΦ	7-34 mm (Free-cut)	97 mm × 105 mm	M	5 pieces
M16	16642	16 mmΦ	12-15 mm	97 mm × 105 mm	M	5 pieces
M19	16643	19 mmΦ	15-18 mm	97 mm × 105 mm	M	5 pieces
M22	16644	22 mmΦ	18-21 mm	97 mm × 105 mm	M	5 pieces
M25	16645	25 mmΦ	21-24 mm	97 mm × 105 mm	M	5 pieces
M28	16646	28 mmΦ	24-27 mm	97 mm × 105 mm	M	5 pieces
L32	16647	32 mmΦ	27-44 mm (Free-cut)	108 mm × 117 mm	L	5 pieces