


Mini Close

Appliances for Bathing

Features

- Recommended when stool volume is less such as after irrigation or during bathing.
- Adhesive ingredient with less skin irritation is used for faceplate.
- Small and oval faceplate reduces uncomfortableness and does not interfere with body motion.
- Flesh-colored bag makes it less noticeable.



Type	Product Code No.	Specification			Contents Per Box
		Starter Hole Size	Max Available Size	Bag Size (LxW)	
M	12641	40 mm Φ	- 39 mm	120 mm x 110 mm	30 pieces
L	12642	40 mm Φ	- 39 mm	140 mm x 130 mm	30 pieces