

HN-0.5 / Breast Report

**Centre Saint-Jean
Saint-Doulchard**

<Challenge>

For breast treatments, we encounter difficulties in repositioning the arms in the arm rest of our breast support (MBLX from Macromedics).


There are many different ways of repositioning the arms inside the arm rests, creating uncertainty about repositioning,

All the more so for breast + lymph node treatments (susclav, subclav, axillary).

<Application>

We use Moldcare HN-0.5 (15 x 20 cm) positioned and molded in the armrests to ensure session after session ideal repositioning of the shoulders and humerus.

This allows us to adapt the size of the armrest to the morphology of the patient's arm.



<Result>

Excellent results.

Moldcare helps us to reposition all our patients on a daily basis, and at the same time to ensure optimized registration, even for highly mobile treatment areas.